

MATERIALI VALIDATI*



MATERIALE	NOME	FABBRICANTE	INTERNET
Registrazione del morso	Metal Bite	r-dental	www.r-dental.com
Resina acrilica autopolimerizzante a freddo <i>Materiale bicomponente</i>	<ul style="list-style-type: none"> ■ ProBase Cold Monomer ■ ProBase Cold Polymer 	Ivoclar	www.ivoclarvivadent.com
Resina radiopaca <i>Materiale bicomponente</i>	<ul style="list-style-type: none"> ■ Acryline X-Ray Powder (25% Ba-SO₄) oppure Acryline X-Ray Powder DVT (10% BaSO₄) ■ Acryline Liquid 	Anaxdent	www.anaxdent.com
Materiale di termoformatura <i>1,5mm – 2mm</i>	Erkodur	Erkodent	www.erkodent.com

*L'utilizzo di materiali non adatti per la realizzazione della dima radiologica può causare rischi per la salute del paziente. Può anche influenzare negativamente le prestazioni e le caratteristiche del prodotto. Utilizzare esclusivamente materiali validati da SICAT.

SICAT
SUPPORT

SGL@SICAT.COM

T +49 228 286 206-620


FABBRICANTE

 SICAT GMBH & CO. KG
 FRIESDORFER STR. 131-135
 53175 BONN
 GERMANIA
 WWW.SICAT.COM

DD30QUG206

VERSIONE:

2024-02-20