

## MATERIALI RACCOMANDATI



MATERIALE	NOME	PRODUTTORE	INTERNET
<b>Resina acrilica autopolimerizzante a freddo</b>	Biocryl Ice	Great Lakes Orthodontic	<a href="http://www.greatlakesortho.com">www.greatlakesortho.com</a>
	ProBase Cold	Ivoclar	<a href="http://www.ivoclarvivadent.com">www.ivoclarvivadent.com</a>
	PalaXpress	Kulzer	<a href="http://www.kulzer.com">www.kulzer.com</a>
<b>Resina acrilica radiopaca</b>	Biocryl X	Great Lakes Orthodontic	<a href="http://www.greatlakesortho.com">www.greatlakesortho.com</a>
	Biocryl XR	Great Lakes Orthodontic	<a href="http://www.greatlakesortho.com">www.greatlakesortho.com</a>
	Acryline X-ray	Anaxdent	<a href="http://www.anaxdent.com">www.anaxdent.com</a>
	Acryline X-ray DVT	Anaxdent	<a href="http://www.anaxdent.com">www.anaxdent.com</a>
<b>Pellicola per termostampaggio da 1,5mm - 2,0mm</b>	Biocryl splint	Great Lakes Orthodontic	<a href="http://www.greatlakesortho.com">www.greatlakesortho.com</a>
	Erkodur	Erkodent	<a href="http://www.erkodent.com">www.erkodent.com</a>
<b>Materiale per la registrazione del morso</b>	Metal Bite	R-Dental	<a href="http://www.r-dental.com">www.r-dental.com</a>
	Flexitime Bite	Kulzer	<a href="http://www.kulzer.com">www.kulzer.com</a>

**SICAT**
**SUPPORT**

SGL@SICAT.COM

T +49 228 286 206-620

**FABBRICANTE**

 SICAT GMBH & CO. KG  
 FRIESDORFER STR. 131-135  
 53175 BONN  
 GERMANIA  
 WWW.SICAT.COM

DD30IFU106

 FECHA:  
 2021-03-11