





SICAT 



START SMART WITH SICAT

Smart 3D solutions – subscription now available!





INNOVATION. VARIETY. SICAT.

About SICAT	4
History	6
Products	8
Implantology	14
Endodontics	18
Sleep Medicine	22
Functional Dentistry	26
Services	30
Contact	31

With products based on 3D X-ray, SICAT paves the ground for advanced digital dentistry.

**MADE IN GERMANY.
MADE IN BONN.
MADE BY SICAT.**



SICAT's expertise covers the development of 3D-based analysis and planning software for various dental disciplines as well as a broad range of surgical guides and therapeutic appliances to directly implement treatment.

More than 25,000 software installations, 250,000 surgical guides and therapeutic appliances so far – for customers worldwide.

SICAT offers solutions, which enhance the digital chairside workflow. Products based on artificial intelligence and optimized for 3D printing underscore SICAT's prowess as an innovation leader in dentistry.



Owner-operated since 2004 – In-house expertise

More than 120 SICAT employees work hand in hand, from development through in-house production, case check to technical support, making SICAT a full service provider for dentists.

High quality standards – MDR certified since 2022

All products meet the requirements for medical devices defined in the laws and regulations.




SICAT is setting trends and introduces world firsts to the market.

SINCE 2002

Development of an optimized clinical workflow based on CBCT


2007

Galileos Implant




2011

Integration of intraoral scans



2014

SICAT Function & SICAT JMT+



GUIDED IMPLANTOLOGY

FULLY DIGITAL GUIDED IMPLANTOLOGY

TMD TREATMENT & 4D RESTORATIONS

2004

Foundation of SICAT GmbH & Co. KG. Since 1.1.2019 100% independent with fully open interfaces.



SICAT CLASSICGUIDE




SICAT OPTIGUIDE
SICAT DIGITALGUIDE



SICAT OPTIMOTION

2016

SICAT Air




**SLEEP APNEA
TREATMENT**



SICAT OPTISLEEP

2018

SICAT Endo



**GUIDED
ENDO**



SICAT ACCESSGUIDE

2021

AI based ready-to-print STL surgical guide design for in-office 3D printing




**ARTIFICIAL INTELLIGENCE &
3D PRINTING**



SICAT *RAPID* DIGITALGUIDE
POWERED BY ARTIFICIAL INTELLIGENCE

Multidisciplinary treatment options with the SICAT Suite to treat patients comprehensively.




TOP 5 REASONS FOR THE SICAT SUITE

- ✓ **Open** for data of all CBCT manufacturers and optical scanner
- ✓ **Simple interface** to CBCT manufacturer software and automatic import of CBCT data into SICAT Suite
- ✓ **Only one-time installation** required, only one-time to import data
- ✓ **Shared data** across all SICAT applications
- ✓ **Comprehensive treatment planning** for various disciplines, optimum patient consultation included

Same workflow for all applications



ANALYSIS



PLANNING



TREATMENT

Seamless from analysis to treatment.

SICAT SMART

*The fully flexible
SICAT Suite
subscription plan.*





***SICAT Smart** allows dental practices of all sizes to economically access SICAT's world of software applications.*



SMART PAYMENT

Using a SICAT Suite application of choice for **only \$99 per month.**



SMART START

Access to the **SICAT eCademy** and all courses of the chosen application(s).



SMART UPDATES

Updates included – always up-to-date with the most current version of SICAT software

START SMART
*with the SICAT
eCademy.*



ecademy.sicat.com





*A **broad range** of video tutorials to get to know the SICAT Suite applications step by step.*

- ✓ **Getting started with SICAT Software** – easier than ever before!
- ✓ **Learning with video tutorials** whenever and from wherever you want. 100% flexible and at your own pace.
- ✓ **Unlimited access** to courses of the chosen application(s) for all SICAT Smart users.




SICATIMPLANT

Implant planning and diagnostics with a comprehensive software solution.

“The SICAT Implant planning software is the most comprehensive and intuitive in the market. Simply put, SICAT Implant’s workflow is logic excellence. I find myself using it even though I can freehand all of my dental implants. The SICAT workflow is clearly spoiling me.”



Dr. Daniel J Vincent DDS




OVER 13,000 implants, abutments and sleeve systems are integrated in the implant database which is always available in the most current version.



SICAT IMPLANT CONVINCES with its modern user interface and ease-of-use.



MORE THAN 20 YEARS OF DEVELOPMENT and experience are reliably reflected in the algorithm when overlaying optical scans with the CBCT data.



Why to use SICAT Implant

3D View

- **Optimized visualization** for precise analysis and diagnosis of anatomical structures
- **Matching** of any CBCT DICOM data with STL, SSI & SIXD
- **Autofocus** for objects in all section planes to easily navigate
- **360° rotation** around focal point for clear positioning and orientation



Implant aligned view

- **Autofocus** of the implant in all planes for easy positioning
- **360° rotation** around implant for detailed view of neighboring anatomical structures
- **Collision warning** between nerve & implant and adjacent implants
- **Outstanding visualization** and rendering of the 3D model for checking the implant planning

*From planning in
SICAT Implant
to treatment with
SICAT surgical
guides.*

“

*From beginner to pro, from stone model
to in-office printing – SICAT surgical
guides are **the** right solution for every
patient situation and practice equip-
ment.*



Dr. Gerhard Zündorf

SICAT SURGICAL GUIDES

The right guide for any patient case.



WHY DENTISTS CHOOSE SICAT SURGICAL GUIDES

- ✓ **Certified accuracy** of less than 500µ for each surgical guide manufactured by SICAT
- ✓ **Analog workflow** for precise implementation of patient cases with severe artefacts
- ✓ **Digital workflow** from central manufacturing by SICAT to AI based design for in-house manufacturing.

SURGICALGUIDES

CLASSICGUIDE

Safe choice
for every patient case

- Also for edentulous jaw
- Certified production accuracy
- With stone model
- Insensitive to metal artefacts

OPTIGUIDE

High quality,
digital workflow

- Tooth-supported
- Certified production accuracy
- Based on optical impressions

DIGITALGUIDE

Case check, support
and design by SICAT

- Turnaround 12 - 48 hours
- Tooth-supported
- For 3D printing in-office or a lab of choice

DIGITALGUIDE

Automatic design powered by
artificial intelligence

- Download available in less than 20 minutes
- Tooth-supported
- For in-office 3D printing

DIGITALGUIDE *POWERED BY ARTIFICIAL INTELLIGENCE*



Download available in
less than 20 minutes



Ready-to-print STL design



Only **\$19.90** for the final
STL design file




Treating minimally-invasive, seeing more.

“When a tooth is to be preserved, I use SICAT Endo to treat successfully and minimally-invasive. Through the new subscription model the financial is manageable. Time savings, greater safety and the patients are impressed.

Dr. Kathrin Haugk

OPEN
for all CBCT
& optical
scanner



THE DETAILED VISUALIZATION of the root canal system in 3D offers the practitioner all-around safety for the upcoming treatment because conventional images often do not show the exact position and number of roots and canals.

IN ADDITION TO THE 3D SCAN, overlaying with 2D images is also an option.

A SURGICAL GUIDE can be ordered directly from within the software.



SUBSCRIPTION
now available



Why dentists benefit from SICAT Endo

Precise evaluation of the patient situation

- Reliable determination of the number of root canals
- Visualization of curved and calcified canals
- Root canal aligned 360° view
- Exact determination of working length and depth of cavity access



Optimizing the chairside workflow

- Less amount of time for chairside treatment thanks to precise evaluation in advance
- Switching between 2D and 3D data with just one click
- Case analysis is transferred from the patient chair to the desktop pc
- Optimized patient communication through visualization and case report


Guided Endo with SICAT – fast, safe & minimal-invasive treatment.

“
*My confidence in the accuracy of the guide
along with my patients' relief is invaluable.
Even when I do not need to make a guide I use
SICAT Endo routinely to map the canal and
compare that measurement to my apex locator
results giving me confidence that I am always
exactly where I want to be in each canal.*



Dr. Scott Resnick

SICAT ACCESSGUIDE & SICAT RAPID ACCESSGUIDE



WITH THE SICAT SURGICAL GUIDES for guided endodontics, root canals can be accessed minimally invasive by both general practitioners and specialists.

EVEN THE TREATMENT of complicated anatomies and obliterations becomes possible.

AI based guide design or ready-to-print guide design

RAPID ACCESSGUIDE

ACCESSGUIDE

\$19⁹⁰



COMING SOON


Same day treatment possible

- Available in less than 25 minutes for 1 tooth
- \$19.90 per tooth, all canals included
- Ready-to-print STL design powered by AI
- Available for orthograde treatment



Accurate, safe and predictable

- Case check & design by SICAT experts
- Safe orthograde access of canals
- Final surgical guide with certified accuracy



Appropriate drill in the SICAT store

- Available in 16mm and 24mm, Ø 1.2mm
- Individual marking of the drill depth by means of a proven rubber stopper
- Compatible sleeves can be ordered in the SICAT Store




Airway analysis, patient communication and treatment planning all-in-one.

“

SICAT Air software and its unique ability to segment the airway and provide a volumetric analysis of the patient's airway is a great visual tool patients could easily relate to and understand. Starting conversation becomes extremely simple, using our routine cone beam scan that was already in the patient's chart.



Dr. Endre Selmeczy



OPEN
for all CBCT
& optical
scanner

SUBSCRIPTION
now available



USING SICAT AIR, the practitioner becomes the key figure for detecting and treating obstructive sleep apnea and snoring.


BASED ON EXISTING CBCT scans, the software provides a quick overview of the patient's airway situation.

THE VISUALIZATION OF THE AIRWAY in treatment position increases the patients' adherence to their SICAT **OPTISLEEP**, which can be ordered immediately.

Key arguments for SICAT Air

Clarifying patient situation

- **Semi-automatic segmentation** of the relevant airway region as well as calculating volume and area of the airway
- **Visualization** of conspicuous structures and colored highlighting
- **Direct analysis** of airway behavior in case of mandibular protrusion
- **Integrated patient education** and PDF handout



Added value for the dental practice

- **New revenue sources** by expanding the patient base to include all those affected by a widespread disease
- **Patient website** www.optisleep.com and flyer for the waiting room



*Successfully
treat snoring
and OSA.*



DIGITAL DATA RECORDING and milling ensure a high accuracy of fit of the appliance.

PATIENTS LOVE the slim design, the good fit as well as the horizontal and vertical freedom of movement.

THANKS TO CONNECTOR SETS in 10 different lengths, the treatment position can be individually adapted to the patient situation.



This is why dentists recommend SICAT OPTISLEEP

- ✓ **High patient comfort** – easy insertion and removal
- ✓ **Long life cycle** – milled from hygienic, odorless, non-staining PMMA
- ✓ **Good tolerance** – biocompatible, free of metals and plasticizers
- ✓ **Short delivery times** – implementation according to individual needs
- ✓ **Identical and cost-effective reproduction** in case of loss – digital data is stored



Fully digital for happy patients.


SICATFUNCTION

Individual functional prosthetics and TMD treatment

“
SICAT Function changed my life 7 years ago! It helped me redefine many dogmas in occlusion, function and hinge axis. For the first time I was able to see that the movement of the joint is a combination between rotation and translation.



Dr. Javier Vasquez



WITH SICAT FUNCTION, the dentist is able to analyze real jaw movement and existing disorders in the lower jaw movement.

THANKS TO ANATOMICALLY ACCURATE 3D visualization combined with real dynamic jaw movement data of the patient, TMD complaints can be prevented or causally treated.

SICAT FUNCTION COMBINES all treatment steps in one software, completely digital and therefore reproducible and traceable any time.



SUBSCRIPTION
now available



Why dentists benefit from SICAT Function

Visualization & analysis of individual jaw movement

- **Digital integration** of three types of data (CBCT data, optical surface scan data and jaw motion data) for treatment planning using real lower jaw movement dynamics instead of average data.
- **Motion traces** for any point of the mandible – acquisition of all relevant data for analysis and treatment planning in one session.



..... Functional prosthetics

- **Functional prosthetics** – patient-individual measuring of all necessary parameters for the programming of virtual articulators.

Patient-individual diagnostics and treatment of temporomandibular dysfunction (TMD).

“The **OPTIMOTION** therapeutic appliance is patient-friendly and is being used by dentists for treating TMD. It is patient-individual and helps to influence the patient's mandibular movements in order to ensure a correct occlusal set-up and regain a correct mandibular posture.



Dr. Manuele Basili



SICAT OPTIMOTION is a therapeutic appliance for the temporary treatment of TMD and the relieve of associated pain.

HIGH ACCURACY OF FIT of the therapeutic appliance due to the completely digital manufacturing process.

THE APPLIANCE DESIGN is based on the basic principles of the Michigan appliance and takes into account the patient-individual, dynamic occlusion.



*This is why dentists appreciate **SICAT OPTIMOTION***

- ✓ **High patient comfort** due to slim design
- ✓ **Multiple options** – maxillary or mandibular, with both canine or anterior canine guidance as well occlusal relief or smooth surface
- ✓ **Long life cycle** – milled from hygienic, odorless, non-staining PMMA
- ✓ **Identical reproduction** in case of loss



Fully digital workflow

SICAT SERVICES

If you choose SICAT, you enjoy extensive support!



Fast solutions in case of technical challenges

SICAT Support reliably assists and quickly implements customer requests.



No more open questions thanks to SICAT trainings

A wide range of training offers turn users into SICAT pros in no time.



Individual case check through certified dental technicians

4 eyes principle to provide safety during patient treatment.

CONTACT

YOUR DIRECT PATH TO SICAT

For general inquiries: +49 (0)228 / 286206-0
SICAT Support: +49 (0)228 / 286206-600
Toll Free (USA & Canada): 1-800-550-9961

Sales: sales@sicat.com
Surgical Guides: sgl@sicat.com
Therapeutic Appliances: support@sicat.com
Software: softwaresupport@sicat.com



SICAT's Managing Directors

“
Our business idea originated in medicine. Our goal was to use 3D imaging techniques to plan surgical procedures in the oral and maxillo-facial area safely and to make them even more predictable with the help of surgical guides. Today SICAT stands for innovation and lives innovation.



Jochen Kusch



Dr. Joachim Hey



SICAT GMBH & CO. KG
FRIESDORFER STR. 131 – 135
53175 BONN
GERMANY

T +49 (0)228 / 286206-0
F +49 (0)228 / 286206-971

INFO@SICAT.COM
WWW.SICAT.COM



207/03-23-EN