



## VALIDATED MATERIALS\*



MATERIAL	NAME	MANUFACTURER	INTERNET
<b>Bite registration material</b>	Metal Bite	r-dental	<a href="http://www.r-dental.com">www.r-dental.com</a>
<b>Cold-curing acrylic 2-component material</b>	<ul style="list-style-type: none"> <li>■ ProBase Cold Monomer</li> <li>■ ProBase Cold Polymer</li> </ul>	Ivoclar	<a href="http://www.ivoclarvivadent.com">www.ivoclarvivadent.com</a>
<b>Radio-opaque acrylic 2-component material</b>	<ul style="list-style-type: none"> <li>■ Acryline X-Ray Powder (25% Ba-SO<sub>4</sub>) or Acryline X-Ray Powder DVT (10% BaSO<sub>4</sub>)</li> <li>■ Acryline Liquid</li> </ul>	Anaxdent	<a href="http://www.anaxdent.com">www.anaxdent.com</a>
<b>Thermoforming sheet 1,5mm – 2mm</b>	Erkodur	Erkodent	<a href="http://www.erkodent.com">www.erkodent.com</a>

\*The use of incorrect materials for the preparation of the radiographic template can lead to health risks for the patient. It can also negatively influence the performance and characteristics of the product. Only use materials that have been validated by SICAT for preparation.

**SICAT** 

**SUPPORT**

SGL@SICAT.COM

T +49 228 286 206-620

 **MANUFACTURER**

SICAT GMBH & CO. KG  
FRIESDORFER STR. 131-135  
53175 BONN  
GERMANY  
[WWW.SICAT.COM](http://WWW.SICAT.COM)

DD30QUG202  
REVISION:  
2024-02-20